



Sterile Disposable Skin Stapler (Zustap)
Sterile Disposable Staple Remover (Zurem)

Mkt By Medzus
MEDZUS Healthcare Pvt. Ltd.
1 Picadilly, Royal Palms, Goregaon (E), Mumbai 65

Customer Care: 022-28735685/ Email id: mail@medzushealth.com/ Website: www.medzushealth.com



**Make it a part
of your team**

Skin Stapler

Sterile Disposable Skin Stapler (Zustap)



Skin Stapler Features:

- ◆ Stapler handle design gives good texture to the product & provides good comfortness to the surgeon's hand.
- ◆ Angled Head design provides clear view to ensure precise staple placement. Also angle head design helps the staple to enter easily in to the tissues.
- ◆ The elaborate design of the staple release mechanism actuates the stapler very efficiently.
- ◆ Easy to use, reduce learning curve.
- ◆ Two types of Skin Stapler suitable for small & big hand with any of four staple sizes i.e. Wide (W), Regular (R), High (H), Large (L) and with any of nine types of staple counts.

Choice! of Stapler:

1.Types:

- ◆ **D-Model:** Suitable for big hands; Male Surgeons.
- ◆ **E-Model:** Suitable for small hands; Female Surgeons.

2.Staple Sizes:

- ◆ **Wide (W):** Wide staple is applicable for general part of body.
- ◆ **Regular (R):** Regular Staple is applicable for thin portion of body. Ex. Face/Head.
- ◆ **High (H):** High staple is applicable for fat portion of body. Ex. Fat person's stomach.
- ◆ **Large (L):** Large staple is applicable for special part of body.

Staple Remover

Sterile Disposable Staple Remover (Zurem)



Staple Remover Features:

- ◆ Special design of the remover's jaw makes the remover comfortable & easy for removing staple.
- ◆ The remover can remove all types of staples.

Choice! of Staple Remover:

Types:

- ◆ T-1 Model: Wide Design
- ◆ T-2 Model: Narrow Design

Technical Specification:

- ◆ Material of staple: 304L Grade stainless steel.
- ◆ Sterilization: Ethylene Oxide gas sterilization.
- ◆ Packaging: Packed with DUPONT Tyvek paper.